

NEVELIA®

BI-LAYER MATRIX

**The choice for aesthetic and functional results
in loss of substances¹**

SKIN EXPERTISE WITH A NEW COLLAGEN TECHNOLOGY

SYMATESE: OUR EXPERTISE IN THE FIELD OF COLLAGEN AND SKIN

SYMATESE GROUP is recognized for the **QUALITY AND EXPERTISE OF ITS TECHNOLOGIES** by some leading companies in cosmetics and dermatology. The key features of NEVELIA® have been developed to achieve dermal reconstruction close to native skin thanks to:

- ✓ Our expertise in the field of **SKIN AND DERMAL REGENERATION**
- ✓ Our **SCIENTIFIC APPROACH** to development and manufacturing
- ✓ Our know-how in **COLLAGEN TRANSFORMATION**

The Collagen extraction and purification are the key for the quality of the implant and guarantee

- Collagen quality
- Product safety
 - Calf hides from animal less than 6 months, coming from safe countries
 - Viral and BSE safety by selection of sourcing and several chemical inactivation processes
 - Conformity to current ISO-22442 and Eur. Pharmacopeia, Monograph N°1482

NEVELIA® Collagen Matrix design and crosslinking is the result of our **expertise in the field of collagen**

- The extraction procedure and the freeze-drying process allow to structure the collagen into a matrix with optimal hydrophily, pore structure and pore size.
- The collagen is then crosslinked to adjust the collagen degradation rate while the dermis is regenerated and therefore optimizes the neodermis quality.


The goal is to obtain an optimal quality of collagen and a matrix structure to allow the neodermis formation.

NEVELIA® BI-LAYER MATRIX

NEVELIA® Bi-Layer Matrix is a sterile medical device consisting of a collagen porous layer to promote and guide regeneration and a reinforced silicone layer acting as a pseudo-epidermis.

This matrix serves as a support for cell infiltration, thus contributes to the natural tissue regeneration process. It is resorbed, becoming a vascularized tissue that is histologically very close to the normal dermis, from 2 to 3 weeks after it is implanted².

The silicone layer is removed after dermal regeneration, at the time of the thin split thickness skin graft.




SYMATESE know-how in collagen extraction and transformation over 30 years has permitted to develop and manufacture successfully collagen-based medical devices and components like bone substitutes, hemostatic compresses or coatings for vascular grafts.

NEVELIA® STRUCTURE

NEVELIA® is a three-dimensional porous matrix of stabilized bovine origin type I Collagen.

NEVELIA® is made of a specific native collagen with a large fibrous proportion to keep cell adhesion signals and mechanical structure to support regeneration.
In vitro tests demonstrate an optimized colonization as fibroblasts recognize collagen fibers.




NEVELIA®
Colonization in the matrix thickness, many cells, strong collagen neosynthesis

OTHER BI-LAYER MATRIX
Colonization mainly on the surface, few cells, poor collagen neosynthesis

Additional NEVELIA® features are:

- ✓ Optimized pore size³ (average pore size 100 µm) and open-cell structure promoting nutrients flow and fibroblasts migration




NEVELIA® MATRIX SEM photography

OTHER BI-LAYER MATRIX SEM photography

- ✓ Crosslinking rate for a balanced absorption / regeneration process
- ✓ No GAG added to keep cells attachment potential

NEVELIA® is composed of medical-grade silicone elastomer reinforced with a polyester material.



NEVELIA®
Reinforced silicone sheeting

OTHER BI-LAYER MATRIX
Non-Reinforced silicone sheeting

NEVELIA® UTILIZATIONS

NEVELIA® is used for dermal regeneration in skin loss, especially in:

- Burns surgery (third and deep second degree burns and burns sequelae)
- Chronic wounds surgery (including leg ulcers and diabetic foot)
- Traumatology
- Skin tumors surgery
- Reconstructive plastic surgery

NEVELIA® can also be used in children.

NEVELIA® is used in combination with a thin split thickness skin graft to recreate skin close to normal skin in terms of function and appearance.

BURNS




CHRONIC WOUNDS



TRAUMATOLOGY



SKIN TUMORS



OTHERS




NEVELIA®, AN ALTERNATIVE TO

- Skin Expansion (Single or sequential)
- Flaps (Local, distant, free ...)
- Split-Thickness or Full-Thickness skin autograft
- Dermal graft (Allografts, Xenografts...)

NEVELIA® DERMAL REGENERATION PROCESS



CONTRACTED SCAR
The surgery can be scheduled as soon as the patient is stabilized.




EXCISION OF THE WOUND
All the dead tissue on and around the graft zone must be removed. For the success of NEVELIA® Bi-layer Matrix, all necrotic and scar tissues must be excised down to the viable tissue.



NEVELIA® BI-LAYER MATRIX APPLICATION
NEVELIA® Bi-Layer Matrix must be cut to fit the excised wound size exactly. The collagen matrix must be in direct contact with the excised wound. The matrix will be fixed in place with surgical staples or sutures.



NEODERM FORMATION
The collagen matrix is quickly colonized by the patient's cells and is gradually replaced by an autologous neoderm.



REMOVAL OF THE SILICONE LAYER
When the neoderm is formed, the silicone layer is removed. The reconstructed dermis has a distinctive orange-yellow or light yellow colour and may present with slightly reddish areas, sign of a good revascularization.



HARVESTING OF THE GRAFT
A thin split thickness skin graft must be taken, if possible from an area similar in colour to the reconstructed area. It must be removed with dermatome and may be meshed.



THIN SPLIT THICKNESS SKIN GRAFT
The graft is placed on the neoderm and it is fixed with staple or suture.



REGENERATED SKIN
Epidermization. Complete healing of the wound.

NEVELIA® BENEFITS

NEVELIA® provides immediate wound closure and promotes dermal regeneration.

NO GAG
Glycosaminoglycan

< ● Good neo-dermis colonization⁴ >

< ● Collagen hemostatic properties
The hemostatic effect may reduce the risk of hematoma formation under or in the device >

NATIVE COLLAGEN
AND ITS 3D OPEN-CELL STRUCTURE


< ● Adequate management & medical care in mobile areas
● Good trophicity especially necessary in areas of friction, pressure, loss of substance, instable scar >

REINFORCED
SILICONE SHEETING

< ● Possible prolonged placement of NEVELIA® for a better management of certain indications especially in chronic wounds or in case of skin graft lack in extended deep burns >


Other bi-layer matrix

NEVELIA®



< ● Reduction of the risk of certain side effects
The matrix will not tear around the staples or sutures
● Neodermis formation over all the surface without any granulation tissue
● Treated area totally covered, thus reducing the risk of infection. >

PRE-HYDRATED MATRIX PLACED BETWEEN 2 RIGID PLASTIC PROTECTIVE SHEETS



1 2 3 4

< ● Ready for use and ease of aseptic handling "no touch" >

8

9

NEVELIA® BI-LAYER MATRIX CLINICAL BENEFITS FOR SPECIFIC INDICATIONS




BURNS SURGERY

- Waiting donor site availability for extended burns




CHRONIC WOUNDS SURGERY

- Good trophicity especially necessary in areas of friction, pressure
- Adapted to prolonged healing time



TRAUMATOLOGY

- Good trophicity especially necessary in areas of friction, pressure
- Adapted to prolonged healing time to combine other treatments



SKIN TUMORS SURGERY

- Wise choice awaiting anatomopathologist results
- Facilitates monitoring of cancer recurrence
- Ease of care for the elderly through faster mobilization

ADDITIONAL BENEFITS

NEVELIA® by-layer matrix can improve treatment in a number of situations.

- ✓ **Immediate availability**
- ✓ **Efficiency** of the surgical technique
- ✓ **Creation of a thicker protective tissue and reestablishment of the gliding plan**
in noble elements coverage like muscles, bones, tendons
- ✓ **Better functional and aesthetic results** compared to common dermo-epidermal graft
- ✓ **No additional scars in traumatology** (versus flaps)
- ✓ **Reduction risk sequelae of donor site**
- ✓ **Does not prevent alternative treatments**
- ✓ **Reduction of hypertrophic scars and keloids occurrence**
- ✓ **Faster procedure** which can be performed in ambulatory room for specific indications (hand, skin tumors)

Children : allows better adaptation of tissue to child growth (less tensions & flanges)

NEVELIA® RESULTS

CLINICAL CASE 1: BURN

- Male, 75 yrs old
- Third-degree burn (left foot)
- Burn date: 07 Jun 2014
- Enrollment date: 01 Jul 2014

Pr Vincent CASOLI
Bordeaux - France



3RD DEGREE BURN



WOUND BED
PREPARATION



NEVELIA® GRAFT



FOLLOW UP



SILICONE REMOVAL,
NEODERMIS FORMED



SPLIT DERMO
EPIDERMAL GRAFT



DAY 73



6 MONTHS




12 MONTHS

CLINICAL CASE 2: CHRONIC WOUND

- Female, 60 yrs old
- Infected chronic wound (right foot)
- Wound date: 28 Nov 2014
- Enrollment date: 17 Feb 2015

Dr Sergiu FLUIERARU
Dr Luc TEOT
Montpellier - France




PREOPERATIVE



WOUND BED
PREPARATION



NEVELIA® GRAFT



SILICONE REMOVAL,
NEODERMIS FORMED



SPLIT DERMO
EPIDERMAL GRAFT



3 MONTHS



6 MONTHS



12 MONTHS

CLINICAL CASE 3: **TRAUMATOLOGY**

- Male, 66 yrs old
- Crush injury (right leg) with loss of soft tissue
- Accident date: 06 Mar 2014
- Enrollment date: 21 May 2014


Pr Vincent CASOLI
Bordeaux - France



WOUND BED
PREPARATION




NEVELIA® GRAFT



FOLLOW UP



SPLIT DERMO
EPIDERMAL GRAFT



DAY 42



6 MONTHS




12 MONTHS

CLINICAL CASE 4: **SKIN TUMORS**

- Female, 62 yrs old
- Wound related to melanoma exeresis (right leg)
- Surgery date: 17 Dec 2014
- Enrollment date: 30 Jan 2014

Pr Vincent CASOLI
Bordeaux - France




WOUND BED
PREPARATION



NEVELIA® GRAFT



FOLLOW UP



SPLIT DERMO
EPIDERMAL GRAFT



6 MONTHS



12 MONTHS

CLINICAL CASE 5: RECONSTRUCTIVE SURGERY

- Male, 57 yrs old
- Open wound on amputation stump (right knee)
- Accident date: 18 Apr 2014
- Enrollment date: 23 May 2014

Pr Vincent CASOLI
Bordeaux - France



BEFORE TREATMENT
PREPARATION OF THE
WOUND WITH NPWT



DAY 0
GRANULATION TISSUE



DAY 0
NEVELIA® GRAFT



DAY 31
FOLLOW UP



DAY 31
SILICONE REMOVAL,
NEODERMIS FORMED



DAY 31
SPLIT DERMO EPIDERMAL
GRAFT



6 MONTHS



12 MONTHS

CLINICAL CASE 6: RECONSTRUCTIVE SURGERY

- Female, 68 yrs old
- Necrotizing dermohypodermatitis (right hand)
- Accident date: 07 Jul 2014
- Enrollment date: 30 Jul 2014

Pr Vincent CASOLI
Bordeaux - France



DAY 0
WOUND BED
PREPARATION



DAY 0
NEVELIA® GRAFT



DAY 12
FOLLOW UP



DAY 22
BEFORE
SILICONE REMOVAL



DAY 22
SILICONE REMOVAL,
NEODERMIS FORMED




DAY 22
SPLIT DERMO EPIDERMAL GRAFT



DAY 22



6 MONTHS



Availability of the device can be subject to local regulatory registration. Please refer to appropriate Instructions for Use for complete product information. This product does not contain latex. Warning: Applicable laws restrict the sale of these devices by or on the order of a physician. NEVELIA® is a trademark of SYMATESE.

NEVELIA® is indicated for dermal regeneration in individuals with skin loss, particularly in the following fields:

- Burns surgery (third and deep second degree burns)
- Reconstructive plastic surgery
- Traumatology

NEVELIA® bi-layer matrix is particularly useful for:

- patients who are unable to supply sufficient donor skin for an autograft at the time of excision
- When the physiological condition of the patient does not allow the autograft

NEVELIA® bi-layer matrix must not be used in patients presenting with:

- clinical signs of wound infection
- an allergic predisposition or known allergy to bovine collagen or silicone

1 Prospective evaluation of NEVELIA® in terms of safety and efficacy for third-degree burns treatment and reconstructive surgery in the University Hospital of Bordeaux: electronic poster presentation at the EWMA 2017 conference in Amsterdam.

2 Instructions for use MCSN01D (May 2013)

3 Yannas, I. V., E. Lee, et al. (1989). «Synthesis and characterization of a model extracellular matrix that induces partial regeneration of adult

mammalian skin.» Proc Natl Acad Sci U S A 86(3): 933-937.

Harley BA, Kim HD, Zaman MH, Yannas IV, Lauffenburger DA, Gibson LJ. Microarchitecture of three-dimensional scaffolds influences cell migration behavior via junction interactions. Biophys J. 2008 Oct;95(8):4013-24.

4 CHAJRA Hanane [18-12-2006] , Mise au point de nouveaux biomatériaux à base de collagène pour la réparation tissulaire cutanée

Development of new collagen-based biomaterials for skin tissue repair



For more information or to place an order,
please contact:



SYMATESE AESTHETICS Customer service
12, rue Jules Verne - 69630 Chaponost - France
Phone: +33(0)4 78 56 72 80
Fax: +33(0)4 78 56 00 48
contact@symatese.com
www.symatese.com